


CONANT
2276.HK





**1.74 high-index
PLUS POWER**

■ Comparison of lens thicknesses of 1.50, 1.60, 1.74 under the same degree

Example: -8.00



/ REAL SHOT



Conant's 1.74 high index hyperopia series uses MR-174™ raw materials from Japan's Mitsui Chemicals, which is a super high index and light lens material. Its excellent optical performance allows users to experience comfort and lightness, at the same time it improves color reproduction and clarity.

■ Lighter and more comfortable



■ Multifunction



■ Derived from plants



■ Luminosity range

Diam	Sph.	Cyl.(-)
65mm	3.00-10.00	-0.25~-4.00
	10.25-12.00	Null

