

Conant®

CONANT
2276.HK



POLARIZED

1.50

1.60

1.67

LENS



HD

Clear view



Real color



Anti glare

UV

UV protection



Comfortable vision



Impact resistant

EXCELLENT PVA COATING

Conant has the most advanced polarized lens production technology. The high quality polarized lenses own an uniform and stable color. The thickness of PVA polarized film is 1/5 of the other ordinary polarized films, and the light transmission effect of our polarized lenses largely exceeds the ordinary polarized lenses' effect. The high temperature softening process makes the film combine with the base material perfectly and is not easy to crack.



COMFORTABLE VISION FEELING

Designed to remove harmful lights and block glare lights, with high definition and high color contrast, Conant's polarized lens gives good appearance and comfortable feelings.



Normal light



Use of Conant lenses

Light Colour Polarized

High transmittance, Restore visual clarity, Visual field is more natural.



Double Combination, All-weather protect eyes.