

Conant[®]

CONANT
2276.HK

**SingleVision
Finished Lenses**





SingleVision Finished Lenses

Finished Lens



1.50

Design.	Diam.(mm)	Sph. (+)	Sph. (-)	Cyl. (-)
Spherical	60	0.25 ~ 8.00	n/a	0.00 ~ 4.00
Spherical	65	0.25 ~ 8.00	0.00 ~ 6.00	0.00 ~ 4.00
Spherical	70	0.25 ~ 4.00	0.00 ~ 6.00	0.00 ~ 4.00
Spherical	75	0.25 ~ 4.00	0.00 ~ 6.00	0.00 ~ 2.00
Spherical	80	n/a	n/a	n/a

1.60

Design.	Diam.(mm)	Sph. (+)	Sph. (-)	Cyl. (-)
Spherical	65	2.25 ~ 8.00	n/a	0.00 ~ 4.00
Spherical	70	0.25 ~ 2.00	6.25 ~ 8.00	0.00 ~ 4.00
Spherical	70	2.25 ~ 8.00	8.25 ~ 10.00	0.00 ~ 2.00
Spherical	75	n/a	0.00 ~ 6.00	0.00 ~ 4.00
Aspheric	65	2.25 ~ 8.00	n/a	0.00 ~ 4.00
Aspheric	70	0.25 ~ 2.00	6.25 ~ 8.00	0.00 ~ 4.00
Aspheric	70	2.25 ~ 6.00	8.25 ~ 10.00	0.00 ~ 2.00
Aspheric	75	n/a	0.00 ~ 6.00	0.00 ~ 4.00

1.67

Design.	Diam.(mm)	Sph. (+)	Sph. (-)	Cyl. (-)
Aspheric	65	2.25 ~ 10.00	n/a	0.00 ~ 4.00
Aspheric	65	n/a	8.25 ~ 10.00	4.25 ~ 6.00
Aspheric	65	n/a	10.25 ~ 12.00	2.25 ~ 4.00
Aspheric	65	n/a	12.25 ~ 14.00	0.25 ~ 2.00
Aspheric	70	n/a	0.00 ~ 8.00	4.25 ~ 6.00
Aspheric	70	0.25 ~ 2.00	6.25 ~ 10.00	0.00 ~ 4.00
Aspheric	70	2.25 ~ 6.00	10.25 ~ 12.00	0.00 ~ 2.00
Aspheric	70	n/a	12.25 ~ 15.00	n/a
Aspheric	75	n/a	0.00 ~ 6.00	0.00 ~ 4.00
Aspheric	75	n/a	6.25 ~ 10.00	0.00 ~ 2.00
Aspheric	80	n/a	0.00 ~ 6.00	0.00 ~ 2.00

1.56

Design.	Diam.(mm)	Sph. (+)	Sph. (-)	Cyl. (-)
Spherical	65	0.25 ~ 6.00	n/a	0.00 ~ 4.00
Spherical	70	0.25 ~ 6.00	8.25 ~ 9.00	0.00 ~ 4.00
Spherical	70	n/a	0.00 ~ 6.00	2.25 ~ 6.00
Spherical	70	n/a	6.25 ~ 7.00	0.00 ~ 6.00
Spherical	70	n/a	7.25 ~ 8.00	0.00 ~ 5.00
Spherical	70	n/a	9.25 ~ 10.00	0.00 ~ 3.00
Spherical	70	n/a	10.25 ~ 11.00	0.00 ~ 2.00
Spherical	70	n/a	11.25 ~ 12.00	0.00 ~ 1.00
Spherical	70	n/a	12.25 ~ 13.00	n/a
Spherical	72	n/a	0.00 ~ 6.00	0.00 ~ 2.00
Aspheric	65	0.25 ~ 6.00	n/a	0.00 ~ 4.00
Aspheric	70	n/a	0.00 ~ 6.00	2.25 ~ 4.00
Aspheric	70	n/a	6.25 ~ 8.00	0.00 ~ 2.00
Aspheric	70	n/a	8.25 ~ 10.00	n/a
Aspheric	72	n/a	0.00 ~ 6.00	0.00 ~ 2.00

1.74

Design.	Diam.(mm)	Sph. (+)	Sph. (-)	Cyl. (-)
Aspheric	65	n/a	10.25 ~ 11.00	5.25 ~ 6.00
Aspheric	65	n/a	11.25 ~ 12.00	4.25 ~ 6.00
Aspheric	65	n/a	12.25 ~ 13.00	3.25 ~ 5.00
Aspheric	65	n/a	13.25 ~ 14.00	2.25 ~ 4.00
Aspheric	65	n/a	14.25 ~ 15.00	1.25 ~ 3.00
Aspheric	65	n/a	15.25 ~ 16.00	0.00 ~ 2.00
Aspheric	65	n/a	16.25 ~ 17.00	0.00 ~ 1.00
Aspheric	65	n/a	17.25 ~ 18.00	n/a
Aspheric	70	n/a	6.25 ~ 10.00	0.00 ~ 6.00
Aspheric	70	n/a	10.25 ~ 11.00	0.00 ~ 5.00
Aspheric	70	n/a	11.25 ~ 12.00	0.00 ~ 4.00
Aspheric	70	n/a	12.25 ~ 13.00	0.00 ~ 3.00
Aspheric	70	n/a	13.25 ~ 14.00	0.00 ~ 2.00
Aspheric	70	n/a	14.25 ~ 15.00	0.00 ~ 1.00
Aspheric	75	n/a	0.00 ~ 6.00	0.00 ~ 6.00